



DigiFix[®]

External Fixator System



SURGICAL TECHNIQUE FOR THE HAND



United States Patents: 8,277,449

8,282,636

9,066,757

FDA 510(k) Approval: K132731

K192465

Surgical Technique Manual

© 2019 Virak Orthopedics LLC

All rights reserved.



Virak Orthopedics LLC
620 Essex Street, Ste 202
Harrison, NJ 07029

(888) 316-6798
VirakOrtho@gmail.com
www.VirakOrtho.com





TABLE OF CONTENTS

Product Information	4
Indications.....	5
Contraindications.....	5
Warnings and Precautions	6
Design Rationale	8
Materials	10
Preoperative Planning.....	10
Fracture-Dislocations	10
Fractures	10
Volar Plate Arthroplasty.....	10
Contractures	11
Joint Arthrodesis	11
Surgical Technique	13
Required Implant and Equipment	14
Optional Implant	14
Axis Pin Insertion.....	14
Distal K-wires Insertion	15
Multiplanar Stabilization	16
Dynamic Mode	16
Static Mode	18
Other Considerations	20
Required Equipment	20
ORDERING INFORMATION	21
IMPLANTS.....	21
INSTRUMENTS.....	21



PRODUCT INFORMATION


Finger injury is common in sports and accidents. Although most of these injuries can adequately be treated by closed reduction, splinting or internal fixation, there are some that are difficult to effectively manage due to bony comminution and/or joint involvement. Post-traumatic arthritis and contracture of the finger may develop if early mobilization and accurate reduction are not achieved. Finger stiffness and contracture may limit usefulness of the hand in grip and grasp, resulting in functional loss.

In an effort to improve the outcome of complex bone and joint injuries of the finger, we developed an



external fixator that attaches to the bony skeleton on the sides of the finger with smooth stainless steel K-wires. In the dynamic mode, the joint can move to allow concentric interphalangeal (IP) range. In the static mode, the DigiFix® allows rigid fixation across the bony segment(s), providing stability to allow healing.

Figure: (Left) Proximal interphalangeal (PIP) joint fracture-dislocation. (Below) PIP fracture-dislocation after percutaneous application of the DigiFix® External Fixator with distraction.



INDICATIONS

The DIGIFIX® External Fixation System is cleared by the FDA to be used in skeletally mature patients in treatment of:

DYNAMIC MODE:

- 1) Complex fracture-dislocations or fracture-subluxation, unstable dislocations, and pilon fractures of the interphalangeal (IP) joint,
- 2) Post-traumatic contracture of the proximal interphalangeal (PIP) joint, and
- 3) Dupuytren's contracture

STATIC MODE:

- 1) Fractures of the phalanges, and
- 2) Interphalangeal (IP) joint arthrodesis.

CONTRAINDICATIONS

- Poor patient compliance.
- Active infection of the digit.
- Severe osteoporosis whereby there is poor purchase of K-wires into the phalanges.



WARNINGS AND PRECAUTIONS

1. Surgeon familiarity with the device, instrumentation, and surgical technique prior to surgery is crucial to proper device installation.
2. K-wire security in bone and device integrity should be routinely checked by the surgeon or hand therapist. Pin tract infections need prompt recognition and treatment, and may require early device removal.
3. As with all percutaneous skeletal fixation, pin care is important in reducing the incidence and severity of pin tract infections.
4. K-wire placement in bone requires accuracy to avoid damage to nerves, blood vessels and tendons. Use caution when handling the sharp tip of the K-wires.
5. Fracture reduction may be compromised and/or the device may be damaged if the patient accidentally hits the hand against an object or catches the device on clothing or bedding. Instruct the patient to use care to protect the hand.
6. For each case, the *Diamond* should be changed only one direction to avoid breakage.
7. The DigiFix® Brackets are designed to be single use. Do not reuse.

POTENTIAL ADVERSE EFFECTS

The following list includes potential complications typically associated with external fixation devices.

- Prolonged healing
- Distraction of the fracture site
- Pin insertion can result in damage to nerves and vessels
- Infection, painful, swollen or inflamed implant site
- Device fracture
- Loosening or dislocation of the implant requiring revision surgery
- Edema
- Loss of range of motion, joint contracture, joint subluxation, and joint dislocation
- Compartment syndrome
- Septic arthritis
- Delayed unions and intractable pain



- Initial condition may persist or recur requiring further treatment
- Replacement of apparatus or components resulting in reoperation
- Pin insertion leading to tissue necrosis
- External components leading to skin pressure
- Allergic reaction(s) to implant material(s)
- Muscle tendon impalement and excessive operative bleeding
- Nonunion/pseudoarthrosis development; persistence or failure of the bone to regenerate satisfactorily
- Loss of bone mass
- Abnormal growth plate development
- Bone fractures of regenerated bone after device removal
- Discrepancy in limb length
- Excessive motion at the fracture site due to improper device set-up
- Heat build-up and bone necrosis with bone sequestration due to rapid drilling of the bony cortex
- Bone deformity
- Thrombosis, late erosion or arteriovenous fistulas
- Osteomyelitis and persistent drainage at wire site after wire removal
- Inability to compress the bone surface due to poorly secured pins seated in the bone

No other risks are known to be associated with the DigiFix® External Fixator.




DESIGN RATIONALE

The unique design of the DigiFix® allows it to function in either the dynamic or static mode. In the dynamic mode, the axis of rotation of the DigiFix® is aligned to the axis of rotation of the joint, and the distal skeletal fixation is collinear to the mid-axis of the bone to maintain concentric joint reduction. Joint distraction (i.e. lengthened state), if needed, is achieved by elongating (crimping) the *Diamond* portion of the device with pliers (**CAUTION: DO NOT USE THE COMPRESSION SPREADER INSTRUMENT TO ELONGATE THE DIAMOND**). The DigiFix® rotates about the axis pin without this pin moving within the bone.


The DigiFix® can also be used in the static mode for phalangeal shaft fractures or joint arthrodesis. Compression (shortened state), if desired, is generated by opening the *Diamond* with the Compression Spreader instrument.

All K-wires, except for one through the *Axis Pin Hole*, can be locked to the DigiFix® Bracket with set screws.



* Used in STATIC application only


Figure: Diagram of the DigiFix® Bracket.



Lengthened State (distraction)

Normal State

Shortened State (compression)



CAUTION: DO NOT USE THE COMPRESSION SPREADER INSTRUMENT TO ELONGATE THE DIAMOND!



MATERIALS

The DigiFix® is manufactured using 6061 T6 Aluminum, which allows radiographic visualization of the bone and joint during insertion and post-operative imaging. The K-wires are Stainless Steel (316 LVM implant grade).

PREOPERATIVE PLANNING

Before placing the DigiFix® on a patient, it is helpful to discuss the rationale and goal of the treatment. Even in the emergency situation, information about the duration of external fixation, potential problems, expected outcome, and complications should be discussed with the patient.

FRACTURE-DISLOCATIONS

For PIP fracture-dislocations, concentric joint reduction and early motion are the key to success. Percutaneous placement of the DigiFix® in the dynamic mode, with distraction, often achieves both these goals without having to open the fracture. Each fracture-dislocation case must be managed individually and use of the DigiFix® does not preclude an open procedure or internal fixation.

FRACTURES

For comminuted shaft fractures of the middle or proximal phalanx, the surgeon should carefully plan fixation placement and determine optimal surgical exposure. Additional screw fixation, tension band wire and/or K-wire fixation can be used in combination with the DigiFix® as needed. Each fracture must be managed individually. Some patients can tolerate early motion, while others may have more swelling and pain precluding early motion.

VOLAR PLATE ARTHROPLASTY

When using the DigiFix® as adjunct to volar plate arthroplasty, the surgeon must still adhere to the Eaton principles. The collateral ligaments should be excised as part of the volar plate arthroplasty procedure. The anchoring suture through the volar plate should not be tied until the DigiFix® is properly placed and secured. Once the DigiFix® is maintaining joint reduction, the suture in the volar plate can be tied over the button on the dorsum of the finger. The PIP joint is held in nearly full extension for the first two weeks. Then, gentle flexion is initiated, incrementally increasing the range of flexion over the subsequent weeks. Try to increase the amount of flexion of the frame by approximately 20°-30° per week. When not moving the joint passively or actively, the joint is held in nearly full extension. It is also important to check lateral X-rays on a weekly basis to ensure reduction of the joint. Before frame removal, try to achieve as much flexion as possible.



CONTRACTURES

For patients with contractures of the joint, the underlying cause should be addressed and the DigiFix® can be used as adjunct treatment. For example, tenolysis may need to be considered as part of the treatment. If there is doubt as to the degree of tendon adherence, it may be advantageous to perform the capsular release first, regaining requisite passive motion, and delay the tenolysis or tendon reconstruction until a later date. There are occasions where a small zone of flexor adherence is present and tenolysis could be performed as a part of the contracture release and external fixator placement.

The extensor mechanism may be attenuated in patients with long-standing flexion contractures of greater than 50°. In this situation, a boutonniere reconstruction may be needed to provide functioning extensor mechanism.

For boutonniere reconstruction, the PIP should be kept in full extension for at least three weeks. From then, a gradual and incremental program of passive flexion is started, achieving full flexion at week six. When not flexing the joint, the fixator should be kept in full extension, especially overnight. Each week, 10-20° of flexion is added to the program. Active flexion and extension should be cautiously initiated to limit attenuation of the reconstructed extensor mechanism. Some permanent extensor lag should be expected in all of these patients.



Elastomer (dental) bands (red arrow) are used for contracture correction. They attach to the *Hooks* on the DigiFix® Brackets at the distal end and to an anchor pin in the proximal phalanx. The bands provide a low load prolong stretch (LLPS) moment about the joint to overcome contractures over time. The bands are placed dorsal to the axis pin when correcting for flexion contracture, and volar to the axis pin for extension contracture (Figure). Duration of the fixator is typically 6-8 weeks.



Elastastomer (dental) Bands
3/8" (or 3/16") 6.5 oz




Figure: For flexion contracture cases, the elastomer bands are placed **dorsal to the axis pin** creating an extension moment about the PIP joint. The load low, prolong stretch pulls the joint into extension over a period of time.




Figure: For extension contracture cases, the elastomer bands are placed **volar to the axis pin** creating a flexion moment about the PIP joint. The load low, prolong stretch pulls the joint into flexion over a period of time.

DUPUYTREN'S CONTRACTURES

Continuous external force, ie low load prolong stretch, on the palmar fascia of Dupuytren cords has been shown to cause the fibroblasts and their cytoskeletal components (collagen fibers, fibrils and microfibrils) to re-align parallel to the stretching force. The DigiFix® with elastomer bands can be utilized to exert the low load prolong extension torque about the PIP joint. The fixator application can be performed in a staged fashion to gain extension prior to the Dupuytren surgery, or as an adjunct at the same time as the choice Dupuytren surgery. The elastomer bands are placed dorsal to the axis pin of the fixator and are replaced weekly. The DigiFix® and bands do not prohibit active or passive finger flexion during treatment. Duration of the fixator is typically 6-8 weeks.

JOINT ARTHRODESIS

In arthrodesis cases, joint preparation is done in the standard fashion prior to placing K-wires for the DigiFix®. The joint is denuded of cartilage, and bone cuts are performed to achieve the desired amount of angulation. The DigiFix® will accommodate up to about 40° of flexion, depending on the individual anatomy.

DigiFix®

STERILE KIT



SURGICAL TECHNIQUE

REQUIRED IMPLANT AND EQUIPMENT

- DigiFix® Sterile Kit
- Wire driver (not included in set)
- Wire cutter (not included in set)
- Parallel Pliers (not included in set)

OPTIONAL IMPLANT

- 0.062" K-wire (not included)

AXIS PIN INSERTION

Using fluoroscopy, a true lateral of the PIP joint is obtained. A 0.054" (or 0.062" for a larger finger) K-wire is inserted through the axis of rotation of the PIP joint (i.e. axis pin). This axis pin should be placed transverse to the long axis of the finger in the coronal plane and parallel to the joint surface. The axis pin should exit the other side of the finger, and the position checked by fluoroscopy on both the lateral and PA views. If the axis pin is poorly placed (i.e. off-axis), the DigiFix® may not function properly in the dynamic mode. **After placement of the axis pin, reduce and provisionally pin the joint, from dorsal to volar (i.e. dorsal block pin).** If necessary, flex the joint to obtain the reduction before pinning. Cut the provisional pin short.




Figure: The axis pin should be placed transverse to the long axis of the finger and parallel to the joint surface at the center of rotation of the PIP joint. On a true lateral view, the axis pin should appear as a dot (white) in the head of P1. After inserting the axis pin, a dorsal block pin is used to keep the joint reduced. A dorsal placement of the blocking pin will keep it out of the way when driving the K-wires into P2 and placing Brackets.

DISTAL K-WIRES INSERTION

Once the axis pin is in place and the joint is reduced, slide the DigiFix® Bracket over the axis pin, through the *Axis Pin Hole*. Next, using the DigiFix® Bracket as a guide, percutaneously insert a 0.045" K-wire through the most proximal of the *Distal K-wire Holes* and drive it across the middle phalanx to exit on the other side of the finger. (For a larger finger, 0.054" K-wire can be used.) This K-wire should be perpendicular to the long axis and in the mid-axial plane of the middle phalanx to ensure that the PIP joint will be concentrically reduced. The phalanges tend to be more dorsal than expected. It is important not to angulate the pin to engage bone, but instead adjust the DigiFix® Bracket to align with the middle phalanx.




Figure: Bracket is slid over the axis pin.





Figure: An 0.045" (or 0.054" for a larger finger) K-wire is inserted through the most proximal of the *Distal K-wire Holes* and driven across the middle phalanx to exit on the opposite side of the finger. Note a dorsal block pin is not shown in these diagrams.

After placing the first K-wire distally in P2, a second 0.045" (same size as the first) K-wire is then inserted in a similar manner through one of the remaining *Distal K-wire Holes*. The hole is chosen such that the second K-wire is in the middle phalanx (not the DIP joint or distal phalanx). This K-wire should also be perpendicular to the long axis and in the mid-axial plane of the middle phalanx. Check to see that the DigiFix® Bracket is not too close to the skin (leave 3-4mm of space). Some swelling is expected to occur in the postoperative period. Use the Hex Driver to tighten the *Distal Set Screw* to lock the K-wires to the DigiFix® Bracket.





Figure: A second K-wire is placed in a similar manner into P2 through one of the remaining *Distal K-wire Holes*.

MULTIPLANAR STABILIZATION

A second DigiFix® Bracket is placed over the K-wires on the other side of the finger for multiplanar (quadrilateral) stabilization. The Bracket is slid over the axis pin and K-wires using the corresponding holes. In a similar manner, the Bracket is kept 3-4mm off the skin to allow for post-operative swelling. The K-wires are cut flush with the outer portion of the DigiFix®, and the *Distal Set Screw* tightened.





Figure: A second DigiFix® Bracket is placed on the opposite side of the finger using the corresponding holes. The distal K-wires are cut flush with the outer portion of the DigiFix® Bracket, and the axis pin is bent and cut short. The K-wires are locked to the Bracket by tightening the *Distal Set Screw*.

DYNAMIC MODE

If not already done, the axis pin is bent 90° over the outside portion of the DigiFix® on each side of the finger to prevent the Bracket from sliding off the pin. If distraction is required, parallel pliers are used to crimp each *Diamond* independently (in the dorsal-volar plane) to achieve the desired amount distraction. In most cases, 1 mm of distraction is sufficient, which is crimping about half-way.



Figure: The Diamond is elongated by crimping with pliers to gain distraction. **DO NOT USE THE COMPRESSION SPREADER INSTRUMENT TO ELONGATE THE DIAMOND!**

Fluoroscopy is used to confirm that the joint is symmetrically distracted on the PA view and concentrically reduced on the lateral view. The PIP joint is passively ranged to ensure that the joint is freely moving and remains reduced throughout the arc of motion.







Figure: Dynamic Mode. Distraction of the joint is achieved by crimping down on the *Diamond* portion of the Bracket. The joint and DigiFix® fixator rotate about the axis pin. Photograph and radiograph in cadaver showing the joint distraction.

STATIC MODE

In certain clinical scenarios, the surgeon may want to have static fixation of the finger. In such cases, with the DigiFix® Brackets already secured on each side of finger, an additional K-wire (can be size 0.045", 0.054" or 0.062") is inserted through the *Proximal K-wire Holes* and bone. The joint or bone segment is held in the desired position during insertion of this K-wire. The K-wire is cut flush with the Bracket and the *Proximal Set Screw* is tightened to lock the K-wire to the Bracket.




Figure: For static mode, an additional K-wire is placed through the *Proximal K-wire Hole* with the bone segments reduced in the desired position. The K-wire is cut flush with the Bracket and the *Proximal Set Screw* is tightened to lock the K-wire to the Bracket.

If compression is required, the *Diamond* is expanded (in the dorsal-volar plane) using the Compression Spreader instrument (available in the instrument tray). A light Providone-Iodine sponge dressing is placed. The dressing may be replaced in 10-14 days, and the hand washed with soap and water. Motion of the **adjacent joints** may be initiated and tailored to the particular needs of the patient and injury pattern.




Figure: Correct use of the Compression Spreader to open the *Diamond* in the dorsal-volar plane to provide compression. Do not compress beyond what is allowed by the instrument.





Figure: Example of a joint arthrodesis. After the DigiFix® is placed in the static mode, the *Diamond* is expanded to provide compression across the fusion site. Photograph and radiograph of a cadaver PIP joint in compression.

DRESSINGS

A light Providone-Iodine sponge dressing is placed between the skin and Brackets, around the K-wires, and the finger wrapped with Coban. The sponge dressing helps with edema control. The dressing may be replaced in 10-14 days, and the hand washed with soap and water. Motion may be tailored to the particular needs of the patient and injury pattern.




Figure: Providone-Iodine sponge dressing cut and placed between the skin and Brackets.



OTHER CONSIDERATIONS

Depending on the actual fracture pattern or bone geometry and size, the DigiFix® Brackets may be reversed in the proximal-distal (*Hook* is proximal) and/or dorsal-volar (*Hook* is volar) directions. Because of the multiplanar (quadrilateral) fixation, **both Brackets must be pointing in the same direction.**

EXTERNAL FIXATOR REMOVAL

REQUIRED EQUIPMENT

- Hex driver (1.5mm)
- Needle holder (or similar tool)
- Wire cutter (if K-wires are bent)

The DigiFix® may be removed in the office or operating room setting by cutting off the bent portion of the K-wire(s), loosening the set screws and sliding the Brackets off the K-wires. Prep the exposed K-wire ends with betadine. The K-wires are removed in a routine fashion with a needle holder or similar tool. Non-restrictive dressing is placed over the K-wire holes in the skin, and the patient is instructed on therapy, as clinically indicated.



ORDERING INFORMATION

IMPLANTS

Catalog #	Description	Minimum Suggested Qty
Sterile Kit:		
1146-STKT	DigiFix® Sterile Kit	1
Non-Sterile Tray:		
1146-3000	DigiFix® Bracket	2
1146-4004	K-wire 1.1mm (0.045")	4
1146-4005	K-wire 1.4mm (0.054")	2

INSTRUMENTS

Catalog #	Description	Tray Qty
1146-0010	Hex driver 1.5mm	1
1146-5000	Compression Spreader	1

CONTACT US:

888-316-6798 (phone/fax)

VirakOrtho@gmail.com (email)

www.VirakOrtho.com

